

Product Datasheet

HERPES SIMPLEX VIRUS 1/2 antibody (orb430158)

Description	Rabbit polyclonal antibody to HERPES SIMPLEX VIRUS 1/2
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Human HSV, strain F.
Target	HERPES SIMPLEX VIRUS 1/2
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	4.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.