

Product Datasheet

HERPES VIRUS 6 B VARIANT antibody (orb430159)

Description	Mouse monoclonal antibody to HERPES VIRUS 6 B VARIANT
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P
Immunogen	Purified HHV6 nucleocapsids.
Target	HERPES VIRUS 6 B VARIANT
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.