

Product Datasheet

HUMAN IMMUNODEFICIENCY VIRUS 1 gp120 antibody (Biotin) (orb430168)

Description	Goat polyclonal antibody against HUMAN IMMUNODEFICIENCY VIRUS 1 gp120 conjugated to Biotin
Species/Host	Goat
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Purified native gp120 from strain IIB.
Target	HUMAN IMMUNODEFICIENCY VIRUS 1 gp120
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Uniprot ID	Q70626
Expiration Date	12 months from date of receipt.