

Product Datasheet

INFLUENZA A H3N2 antibody (orb430177)

Description	Goat polyclonal antibody to INFLUENZA A H3N2
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Influenza A, Texas 1/77 strain (H3N2).
Target	INFLUENZA A H3N2
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.