

Product Datasheet

INFLUENZA B antibody (orb430181)

Description	Goat polyclonal antibody to INFLUENZA B
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Yamagata strain of Influenza B
Target	INFLUENZA B
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	4.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.