

Product Datasheet

OCHRATOXIN A antibody (orb430290)

Description	Rabbit polyclonal antibody to OCHRATOXIN A
Species/Host	Rabbit
Reactivity	Fungi
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Ochratoxin A conjugated to BSA
Target	OCHRATOXIN A
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.