

Product Datasheet

RABIES VIRUS antibody (orb430372)

Description	Mouse monoclonal antibody to RABIES VIRUS
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	Purified virions of SAD-Vnukovo strain followed by Pitman-Moore strain.
Target	RABIES VIRUS
Preservatives	Ig fraction - liquid
Form/Appearance	Ig fraction - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.