

Product Datasheet

RESPIRATORY SYNCYTIAL VIRUS antibody (orb430374)

Description	Goat polyclonal antibody to RESPIRATORY SYNCYTIAL VIRUS
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FA, ICC, IF, IHC-Fr
Immunogen	Human RSV isolate.
Target	RESPIRATORY SYNCYTIAL VIRUS
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	4.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.