

Product Datasheet

RESPIRATORY SYNCYTIAL VIRUS antibody (Biotin) (orb430375)

Description	Goat polyclonal antibody against RESPIRATORY SYNCYTIAL VIRUS conjugated to Biotin
Species/Host	Goat
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, FA, IHC-Fr
Immunogen	Human RSV isolate.
Target	RESPIRATORY SYNCYTIAL VIRUS
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	4.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.