

Product Datasheet

ROTAVIRUS antibody (orb430382)

Description	Sheep polyclonal antibody to ROTAVIRUS
Species/Host	Sheep
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF
Immunogen	Five serotypes of group A Rotavirus strains Wa (serotype I), DS-1 (serotype II), Ito (III), HOCHI (IV) and 69M (VIII).
Target	ROTAVIRUS
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	5.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.