

Product Datasheet

VACCINIA VIRUS antibody (orb430477)

Description	Rabbit polyclonal antibody to VACCINIA VIRUS
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IH-Fr
Immunogen	Vaccinia virus, New York City Board of Health (NYCBOH) strain.
Target	VACCINIA VIRUS
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	4 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.