

## Product Datasheet

### EXODUS-2 antibody (Biotin) (orb430508)

<b>Description</b>	Rabbit polyclonal antibody against EXODUS-2 conjugated to Biotin
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Recombinant mouse Exodus-2.
<b>Target</b>	EXODUS-2
<b>Preservatives</b>	Purified IgG conjugated to Biotin - lyophilized
<b>Form/Appearance</b>	Purified IgG conjugated to Biotin - lyophilized
<b>Concentration</b>	0.1mg/ml
<b>Storage</b>	-20°C only
<b>Note</b>	For research use only
<b>Isotype</b>	Polyclonal IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P84444</a> , <a href="#">P86792</a> , <a href="#">P86793</a>
<b>Expiration Date</b>	12 months from date of receipt.