

Product Datasheet

MIP-2 antibody (Biotin) (orb430540)

Description	Rabbit polyclonal antibody against MIP-2 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, IH-Fr, IHC-P, WB
Target	MIP-2
Preservatives	Purified IgG conjugated to Biotin - lyophilized
Form/Appearance	Purified IgG conjugated to Biotin - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Uniprot ID	P10889
Expiration Date	12 months from date of receipt.