

Product Datasheet

TRIAD3A antibody (orb430613)

Description	Rabbit polyclonal antibody to TRIAD3A
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A 15 amino acid peptide located near mouse TRIAD3A amino-terminus.
Target	TRIAD3A
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P58283
Expiration Date	12 months from date of receipt.