

Product Datasheet

DYKDDDDK TAG antibody (orb430674)

Description	Rabbit polyclonal antibody to DYKDDDDK TAG
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Keyhole limpet haemocyanin conjugated Enterokinase cleavage site peptide DYKDDDDK. Residues of cysteine and glycine were used to facilitate coupling at the C-terminal end.
Target	DYKDDDDK TAG
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.