

Product Datasheet

HA-TAG antibody (orb430675)

Description	Rabbit polyclonal antibody to HA-TAG
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Keyhole limpet hemocyanin conjugated epitope tag peptide YPYDVPDYA (114-122) from haemagglutinin influenza. A cysteine residue was used to facilitate coupling at the C-terminal end.
Target	HA-TAG
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q03909
Expiration Date	12 months from date of receipt.