

Product Datasheet

VSV-G TAG antibody (orb430773)

Description	Rabbit polyclonal antibody to VSV-G TAG
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide corresponding to the carboxy terminal of vesicular stomatitis virus glycoprotein (VSV-G)
Target	VSV-G TAG
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.