

Product Datasheet

VSV-G TAG antibody (Biotin) (orb430774)

Description	Rabbit polyclonal antibody against VSV-G TAG conjugated to Biotin
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, ICC, IF
Immunogen	Synthetic peptide corresponding to the carboxy terminal of vesicular stomatitis virus glycoprotein (VSV-G)
Target	VSV-G TAG
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.