

## Product Datasheet

### NFkB p65 (phospho-Ser536) antibody (orb430851)

<b>Description</b>	Rabbit polyclonal antibody to NFkB p65 (phospho-Ser536)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-P, WB
<b>Immunogen</b>	Synthetic peptide corresponding to residues surrounding phosphoserine 536 of human NFkB p65 protein.
<b>Target</b>	NFkB p65 (pSer536)
<b>Preservatives</b>	Purified IgG - liquid
<b>Form/Appearance</b>	Purified IgG - liquid
<b>Concentration</b>	0.5mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Polyclonal IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<a href="#">Q04206</a>
<b>Expiration Date</b>	12 months from date of receipt.