

Product Datasheet

WNV antibody (orb430981)

Description	Rabbit polyclonal antibody to WNV
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Immunogen	A peptide corresponding to an 15 amino acid sequence from near the carboxy terminus of West Nile Virus core protein.
Target	WNV
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Core protein (C-Terminal)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P06935
Expiration Date	12 months from date of receipt.