

Product Datasheet

HUMAN IMMUNODEFICIENCY VIRUS 1 gp120 antibody (orb431204)

Description	Goat polyclonal antibody to HUMAN IMMUNODEFICIENCY VIRUS 1 gp120
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant gp160 comprising the C terminus of gp120 and most of gp41.
Target	HUMAN IMMUNODEFICIENCY VIRUS 1 gp120
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	Q70626
Expiration Date	12 months from date of receipt.