

Product Datasheet

HUMAN IMMUNODEFICIENCY VIRUS 1 gp41 antibody (orb431207)

Description	Goat polyclonal antibody to HUMAN IMMUNODEFICIENCY VIRUS 1 gp41
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide CH3CO-YTSLIHSLEESQNQQEKNEQELLELDKWASLWNWF-NH2
Target	HUMAN IMMUNODEFICIENCY VIRUS 1 gp41
Preservatives	Purified IgG - liquid
Form/Appearance	Phosphate buffered saline, 0.1% Sodium Azide (NaN3)
Concentration	p>1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	> 98%
Uniprot ID	Q70626
Dilution Range	ELISA: 1:500-10,000, WB: 1:500-10,000
Expiration Date	12 months from date of receipt.