

Product Datasheet

CCL3L1 antibody (Biotin) (orb431365)

Description	Rabbit polyclonal antibody against CCL3L1 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant pig CCL3L1
Target	CCL3L1
Preservatives	Purified IgG conjugated to Biotin - liquid
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Uniprot ID	F1S1A1
Expiration Date	12 months from date of receipt.