

Product Datasheet

BRCA1 (phospho-Ser988) antibody (orb431546)

Description	Rabbit polyclonal antibody to BRCA1 (phospho-Ser988)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding serine 988 of human BRCA1
Target	BRCA1 (pSer988)
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P38398
Expiration Date	12 months from date of receipt.