

Product Datasheet

GSK3 ALPHA (phospho-Ser21) antibody (orb431607)

Description	Rabbit polyclonal antibody to GSK3 ALPHA (phospho-Ser21)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding serine 21 of human GSK3 ALPHA
Target	GSK3 ALPHA (pSer21)
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P49840
Expiration Date	12 months from date of receipt.