

Product Datasheet

GSK3 BETA (phospho-Tyr216) antibody (orb431609)

Description	Rabbit polyclonal antibody to GSK3 BETA (phospho-Tyr216)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding tyrosine 216 of human GSK3 beta
Target	GSK3 BETA (pTyr216)
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P49841
Expiration Date	12 months from date of receipt.