

## Product Datasheet

### LCK (phospho-Tyr394) antibody (orb431622)

<b>Description</b>	Rabbit polyclonal antibody to LCK (phospho-Tyr394)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Phospho specific-peptide corresponding to residues surrounding tyrosine 394 of human LCK
<b>Target</b>	LCK (pTyr394)
<b>Preservatives</b>	Purified IgG - liquid
<b>Form/Appearance</b>	Purified IgG - liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Polyclonal IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<a href="#">P06239</a>
<b>Expiration Date</b>	12 months from date of receipt.