

Product Datasheet

CONNEXIN 43 (phospho-Ser368) antibody (orb432031)

Description	Rabbit polyclonal antibody to CONNEXIN 43 (phospho-Ser368)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide corresponding to an amino acid sequence within connexin43 which includes phosphorylated Ser368
Target	CONNEXIN 43 (pSer368)
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Uniprot ID	P08050
Expiration Date	12 months from date of receipt.