

Product Datasheet

CD161 / NK1.1 antibody (FITC) (orb432327)

Description	Mouse monoclonal antibody against CD161 / NK1.1 conjugated to FITC
Species/Host	Mouse
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Spleen and bone marrow cells from CE mice.
Target	CD161 / NK1.1
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P27812 , P27814
Expiration Date	12 months from date of receipt.