

Product Datasheet

V5-TAG antibody (orb432439)

Description	Mouse monoclonal antibody to V5-TAG
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, ICC, IF, IHC-Fr, IP, WB
Immunogen	Paramyxovirus Simian-Virus 5 (SV5)
Target	V5-TAG
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P11207
Expiration Date	12 months from date of receipt.