

Product Datasheet

HISTIDINE TAG antibody (orb432482)

Description	Mouse monoclonal antibody to HISTIDINE TAG
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	PAX6 transcription factor linked to histidine tag.
Target	HISTIDINE TAG
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.