

## Product Datasheet

### HISTIDINE TAG antibody (HRP) (orb432483)

<b>Description</b>	Mouse monoclonal antibody against HISTIDINE TAG conjugated to HRP
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	PAX6 transcription factor linked to histidine tag.
<b>Target</b>	HISTIDINE TAG
<b>Preservatives</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.