

Product Datasheet

CALICIVIRUS antibody (Biotin) (orb432962)

Description	Mouse monoclonal antibody against CALICIVIRUS conjugated to Biotin
Species/Host	Mouse
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Target	CALICIVIRUS
Preservatives	Purified IgG conjugated to biotin - liquid
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.