

Product Datasheet

CD200 antibody (FITC) (orb433035)

Description	Rat monoclonal antibody against CD200 conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Mouse CD200-rat CD4 fusion protein.
Target	CD200
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	O54901
Expiration Date	12 months from date of receipt.