

Product Datasheet

CD177 antibody (FITC) (orb433121)

Description	Mouse monoclonal antibody against CD177 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Human granulocytes.
Target	CD177
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid.
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid.
Concentration	0.1mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q8N6Q3
Expiration Date	12 months from date of receipt.