

## Product Datasheet

### CD16/CD56 antibody (FITC) (orb433161)

<b>Description</b>	Mouse monoclonal antibody against CD16/CD56 conjugated to FITC
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC
<b>Target</b>	CD16/CD56
<b>Preservatives</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Form/Appearance</b>	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Concentration</b>	0.1mg/ml
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Isotype</b>	Cocktail
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P08637</a> , <a href="#">O75015</a> , <a href="#">P13591</a>
<b>Expiration Date</b>	12 months from date of receipt.