

## Product Datasheet

### Granzyme B antibody (Biotin) (orb433220)

<b>Description</b>	Mouse monoclonal antibody against Granzyme B conjugated to Biotin
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Purified human Granzyme B.
<b>Target</b>	GRANZYME B
<b>Preservatives</b>	Purified IgG - liquid
<b>Form/Appearance</b>	Purified IgG - liquid
<b>Concentration</b>	0.5 mg/ml
<b>Storage</b>	+4°C, -20°C if preferred
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P10144</b>
<b>Expiration Date</b>	12 months from date of receipt.