

## Product Datasheet

### CD49b antibody (FITC) (orb433230)

<b>Description</b>	Rat monoclonal antibody against CD49b conjugated to FITC
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC
<b>Immunogen</b>	RIL-2 propagated NK1.1+ cells from C57BL/6 mice.
<b>Target</b>	CD49b
<b>Preservatives</b>	Purified IgM conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Form/Appearance</b>	Purified IgM conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
<b>Concentration</b>	0.1 mg/ml
<b>Storage</b>	+4°C DO NOT FREEZE
<b>Note</b>	For research use only
<b>Isotype</b>	IgM
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q62469</a>
<b>Expiration Date</b>	12 months from date of receipt.