

Product Datasheet

CD146 antibody (FITC) (orb433240)

Description	Mouse monoclonal antibody against CD146 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Recombinant human MUC18 (D1-D5) Fc protein.
Target	CD146
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P43121
Expiration Date	12 months from date of receipt.