

Product Datasheet

MOUSE ANTI SHEEP CD45:RPE (orb433348)

Description	Mouse monoclonal antibody against CD45 conjugated to RPE
Species/Host	Mouse
Reactivity	Sheep
Conjugation	RPE
Tested Applications	FC
Immunogen	Ovine efferent lymphatic duct lymphocytes.
Target	CD45
Preservatives	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Clone Number	1.11.32
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)