

Product Datasheet

CD206 antibody (FITC) (orb433362)

Description	Rat monoclonal antibody against CD206 conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Chimaeric CRD4-7-Fc protein
Target	CD206
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q61830
Expiration Date	12 months from date of receipt.