

Product Datasheet

CD41 antibody (Biotin) (orb433395)

Description	Rat monoclonal antibody against CD41 conjugated to Biotin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	Purified murine platelets
Target	CD41
Preservatives	Purified IgG conjugated to biotin - liquid
Form/Appearance	Purified IgG conjugated to biotin - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9QUM0
Expiration Date	12 months from date of receipt.