



Product Datasheet

BRCA1 (Phospho-Ser988) Antibody (orb43343)

Description	BRCA1 (Phospho-Ser988) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P38398
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of H...



Western blot analysis of BRCA1(Phospho-S...