


Product Datasheet

Cyclin E1 (phospho-Thr395) Antibody (orb43349)

| | |
|----------------------------|---|
| Description | Cyclin E1 (phospho-Thr395) Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P24864 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of Cyclin E1(phosp...