


Product Datasheet

Cdc25B (Phospho-Ser149) Antibody (orb43351)

Description	Cdc25B (Phospho-Ser149) Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P30306
Expiration Date	12 months from date of receipt.



Western blot analysis of Cdc25B(Phospho-
...
...