

Product Datasheet

CD172a antibody (FITC) (orb433528)

Description	Mouse monoclonal antibody against CD172a conjugated to FITC
Species/Host	Mouse
Reactivity	Porcine
Conjugation	FITC
Tested Applications	FC
Immunogen	Porcine alveolar macrophages.
Target	CD172a
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q5K4Q3
Expiration Date	12 months from date of receipt.