


Product Datasheet

CK2α (phospho-Thr360/Ser362) Antibody (orb43356)

Description	CK2α (phospho-Thr360/Ser362) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P68400
Expiration Date	12 months from date of receipt.



Western blot analysis of CK2α(Phospho-Thr360/Ser362) antibody.