

Product Datasheet

EPSTEIN-BARR VIRUS LMP2A antibody (orb433762)

Description	Rat monoclonal antibody to EPSTEIN-BARR VIRUS LMP2A
Species/Host	Rat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Bacterial TrpE-LMP2A fusion protein.
Target	EPSTEIN-BARR VIRUS LMP2A
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.