

Product Datasheet

FELINE HERPES VIRUS TYPE 1 antibody (orb433800)

Description	Mouse monoclonal antibody to FELINE HERPES VIRUS TYPE 1
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FELINE HERPES VIRUS TYPE 1
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Expiration Date	12 months from date of receipt.