

Product Datasheet

FELINE LEUKAEMIA VIRUS p15E antibody (orb433880)

Description	Mouse monoclonal antibody to FELINE LEUKAEMIA VIRUS p15E
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Target	FELINE LEUKAEMIA VIRUS p15E
Preservatives	Purified IgG - liquid
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q85525
Expiration Date	12 months from date of receipt.